INNOVATION: INSOLE SYSTEM WITH UNIQUE ANATOMY LAST!

New insole technology relieves the whole Plantar fascia - for a single - and multi combination care

EsConclusion®

indications: - heel spur - Plantar fascia rupture - Arthrosis metatarsophalangeal joint

P P P

Innovative Relief last for the complete inflammation field


Pressure release of the Plantar fascia total

soft bedding in the heel area

Qualität made in German

ESCO Orthopädie Service GmbH Ihre Einlagenmanufaktur

The treatment of heel spur with prefab orthotics gets handled traditionally with relieving pads in the area of the spur. The patient receives a decompression during walk and will be freee of pain. Because often the inflammation, caused by the spur under the calcaneus, a Plantarfascia rupture goes along – and the inflammation field eradiates onto the longitudinal arch of the foot, there is a massive easing, that is based on the correct, anatomical bedding of the Plantar fascia.



The new insole series offers two more professional cares, that redistribute the loading pressure through an anatomic bedding of the Plantar fascia the arch of strain from sinew.

This has the effect of a precise better metabolism care and correspondingly a pain release. This combination is a spupport for a prompt process to free from pains for the patient. Surveys prove it: After 2–3 weeks, the acute symptoms were corrected and the quality and happiness the movement returned. A multifunction ist guranteed through this kind oft treatment because the typical foot deformities can be served optimal by the most existing symptons, which the spur raises. Our innovation ist the technical implementation of the care philosophy into the Carbon insoles technology.

Furthermore to the named indications, the arthritis of the joints from metatarsophalangeal joint can be medicated very well. The known EsCo RigiTec – insole technology is supplemented technically by Plantar fascia bedding.

The foot gets reinforced in its joints. An immediate pain reduction is detectable through the move reduction.

To offer and to realize a variety of therapy-concepts for the orthopaedic technician, we have defined different base materials. These materials are available for the orthopaedic specialist for a choice of his treatment strategy and aditionally an ESD-treatment inside the work safety shoe.





EsConclusion FSP | Carbon

narrow bedding, medium-width, full-length, cover Clarino red on Nora-Lunasoft white

art.–no. sizes 119422-10004 35-48

EsConclusion FSP II Carbon

wide bedding, medium width, full-length, cover Clarino green on EsCo SoftFlex green

art.–no. sizes 119522-10001 35-48



EsConclusion FSP III Carbon

narrow bedding, medium-width, full-length, cover Clarino red on Nora-Lunasoft white

art.–no. sizes 119422-10005 35-48




EsConclusion FSP IV Carbon

wide bedding, medium width, full-length, cover Clarino green on EsCo SoftFlex green

art.–no. sizes 119522-10002 35-48





EsConclusion FSP | HighFlex

narrow bedding, medium width, full-length, cover Clarino black on Dynoform red

artno.	415222-10000
sizes	35-48


EsConclusion FSP II HighFlex

wide bedding, medium width, full-length, cover Clarino blue on Dynoform red

art.-no. sizes

415322-10000 35-48





EsConclusion FSP | EsCombine

narrow bedding, medium width, full-length, cover Clarino black on Dynoform red

art.-no. sizes

124222-10000 35-48

EsConclusion FSP II EsCombine

wide bedding, medium width, full-length, cover Clarino blue on Dynoform red

art.-no. sizes

124322-10000 35-48



EsConclusion FSP | Soft

narrow bedding, medium width, full-length, cover Clarino black on Dynoform red



EsConclusion FSP II Soft

wide bedding, medium width, full-length, cover Clarino blue on Dynoform red

art.-no. sizes 361222-10000 35-48 art.–no. sizes 361322-10000 35-48

EsConclusion HighFlex



EsConclusion HighFlex

medium width, full-length, with cushioned heelspur and Plantar fascia recess, flat



EsConclusion HighFlex

medium width, full-length, with lowered cushioned heel spur- and Plantar fascia recess, cover onSteam grren on EsCo SoftFlex yellow

art.–no. sizes 414022-11229 35-48 art.–no. sizes 414022-10673 35-48

EsConclusion HighFlex





EsConclusion HighFlex

medium width, full-length, with heel spur und Plantar fascia recess, flat, cover onSteam orange on EsCo Soft Flex yellow, metatarsal pad

EsConclusion HighFlex

medium width, full-length, with cushioned heel spur and Plantar fascia recess, raised, cover onSteam anthracite perforated on EsCo SoftFlex yellow, metatarsal pad

art.–no. sizes 414022-10480 35-48

art.–no. sizes 414022-10305 35-48



EsConclusion HighFlex

medium width, full-length, with cushioned heel spur recess, cover onSteam Stone perforated on EsCo SoftFlex yellow, metatarsal pad

art.–no. sizes 414022-10455 35-48